

Dexcom

**Time in Range: An Intuitive CGM
Metric to Integrate in Your
Clinical Practice**

Presented by: Karen Earle, MD



Speaker:
Karen Earle, MD
Chief of Endocrinology and Director of Diabetes Services,
California Pacific Medical Center,
Chief of Division of Medical and Surgical Specialties,
Sutter West Bay Medical Group




Moderator:
Deborah Greenwood, PhD, RN, BC-ADM, CDCES, FADCES
Clinical Education Team
Dexcom

Polling question 1

- What is your current experience with CGM?
 - A. I don't know anything yet, but I saw the Nick Jonas commercial
 - B. I know the basics of CGM but I don't use with patients on a regular basis
 - C. I've worked with CGM several times and feel pretty comfortable
 - D. Expert – I use CGM often and I am able to interpret AGP reports

Current Model of Care






International Consensus on Time in Range Report

A Team Effort

- American Diabetes Association
- American Association of Clinical Endocrinologists
- Association of Diabetes Care and Education Specialists
- European Association for the Study of Diabetes
- Foundation of European Nurses in Diabetes
- International Society for Pediatric and Adolescent Diabetes
- Juvenile Diabetes Research Foundation
- Pediatric Endocrine Society

Diabetes Care Volume 42, August 2019 1593



Clinical Targets for Continuous Glucose Monitoring Data Interpretation: Recommendations From the International Consensus on Time in Range

Diabetes Care 2019;42:1593–1603 | <https://doi.org/10.2337/dci19-0028>

Tadej Battelino,¹ Thomas Danne,² Richard M. Bergenstal,³ Stephanie A. Amiel,⁴ Roy Beck,⁵ Torben Biester,² Emanuele Bosi,⁶ Bruce A. Buckingham,⁷ William T. Cefalu,⁸ Kelly L. Close,⁹ Claudio Cobelli,¹⁰ Eyal Dassau,¹¹ J. Hans DeVries,^{12,13} Kim C. Donaghue,¹⁴ Klemen Dovc,¹ Francis J. Doyle III,¹¹ Satish Garg,¹⁵ George Grunberger,¹⁶ Simon Heller,¹⁷ Lutz Heinemann,¹⁸ Irl B. Hirsch,¹⁹ Roman Hovorka,²⁰ Weiping Jia,²¹ Olga Kordonouri,² Boris Kovatchev,²² Aaron Kowalski,²³ Lori Laffel,²⁴ Brian Levine,⁹ Alexander Mayarov,²⁵ Chantal Mathieu,²⁶ Helen R. Murphy,²⁷ Revital Nimri,²⁸ Kirsten Nrgaard,²⁹ Christopher G. Parkin,³⁰ Eric Renard,³¹ David Rodbard,³² Banshi Saboo,³³ Desmond Schatz,³⁴ Keaton Stoner,³⁵ Tatsuiiko Urakami,³⁶ Stuart A. Weinzimer,³⁷ and Moshe Phillip^{28,38}

Improvements in sensor accuracy, greater convenience and ease of use, and expanding reimbursement have led to growing adoption of continuous glucose monitoring (CGM). However, successful utilization of CGM technology in routine clinical practice remains relatively low. This may be due in part to the lack of clear and agreed-upon glycemic targets that both diabetes teams and people with diabetes can work toward. Although unified recommendations for use of key CGM metrics have been established in three separate peer-reviewed articles, formal adoption by diabetes professional organizations and guidance in the practical application of these metrics in clinical practice have been lacking. In February 2019, the Advanced Technologies & Treatments for Diabetes (ATTD) Congress convened an international panel of physicians, researchers, and individuals with diabetes who are expert in CGM technologies to address this issue. This article summarizes the ATTD consensus recommendations for relevant aspects of CGM data utilization and reporting among the various diabetes populations.

This international consensus report has been endorsed by the American Diabetes Association, American Association of Clinical Endocrinologists, American Association of Diabetes Educators, European Association for the Study of Diabetes, Foundation of European Nurses in Diabetes, International Society for Pediatric and Adolescent Diabetes, JDRF, and Pediatric Endocrine Society.

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)



May underestimate or overestimate an individual's average glucose
(example: A1C of 7% could represent a range between 123-185 mg/dL)

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

2

Does not indicate the extent or timing of hypoglycemia or hyperglycemia

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

3

Does not reveal glycemic variability

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

4

Limited utility for insulin dosing decisions

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

5

Unreliable in patients with hemolytic anemia, hemoglobinopathies, or iron deficiency

Need For Metrics Beyond A1C

A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

6

Underestimates in those with end stage kidney disease or during pregnancy

Need For Metrics Beyond A1C



A1C, %	mg/dL	95% CI
5	97	(76-120)
6	126	(100-152)
7	154	(123-185)
8	183	(147-217)
9	212	(170-249)
10	240	(193-282)
11	269	(217-314)
12	298	(240-347)

7



Correlation with mean glucose can vary among races

Same A1C, but CGM Patterns Drive Different Treatment Plans

A1C 7%




A1C 7%



Adapted from <https://diatribe.org/BeyondA1c>, Assessed March 18, 2021

Continuous Glucose Monitoring Shows Variability



Polling question 2

- What is the % TIR (time in range) goal and % TBR (Time below range) goal for most people with T1D and T2D
 - A. > 70% TIR and < 8% TBR
 - B. > 80% TIR and < 4% TBR
 - C. > 65% TIR and < 5% TBR
 - D. > 70% TIR and < 4% TBR

Core CGM Metrics and Goals for Time in Range (TIR)^{1,2}

Key Metrics

Number of Days with CGM Data

14+ days recommended

Percentage of Time CGM is Active

>70% of data recommended

Mean Glucose

The average glucose


Glucose Management Indicator (GMI)

Approximate A1C levels based on average glucose measured using CGM values


Coefficient of Variation (CV)

Measure of glycemic variability (standard deviation/mean) ≤36% is recommended

T1D & T2D



Older/High Risk T1D & T2D





1. Battelino T et al. *Diabetes Care*. 2019;42(8):1593-1603. 2. American Diabetes Association. *Diabetes Care* 2021;44(Suppl. 1):S73-S84 | <https://doi.org/10.2337/dc21-S006>.

*Includes percentage of values >250 mg/dL **Includes percentage of values <54 mg/dL

Time in Range Targets and Hours per Day^{1,2}

Recommended Time in Range for most people with T1D & T2D



T1D=type 1 diabetes, T2D=type 2 diabetes

1. Battelino T et al. *Diabetes Care*. 2019;42(8):1593-1603. . 2. American Diabetes Association. *Diabetes Care* 2021;44(Suppl. 1):S73-S84 | <https://doi.org/10.2337/dc21-S006>.
 2. *Includes percentage of values >250 mg/dL **Includes percentage of values <54 mg/dL

GOAL



**Increase TIR while
Decreasing TBR**

See **more Green** and **Less Red** (MGLR) on the TIR bar

Glycemic Variability


GLYCEMIC VARIABILITY is the amplitude, frequency and duration of glucose fluctuations¹

Measured by the Coefficient of Variation¹ (CV) Standard deviation/mean=CV

Targets: CV < 36% for most people¹

<33% for people at high risk for hypoglycemia (insulin or sulfonylureas)²

Predictor of hypoglycemia¹



Polling Question 3

- What is the % increase in Time in Range (TIR) that is considered clinically meaningful?
 - A. 1%
 - B. 3%
 - C. 5%
 - D. 10%

Estimated A1C for a TIR level

Each **5% increase** in TIR is clinically significant.¹

Vigersky and McMahon²

(n = 1,137 participants with type 1 and type 2 diabetes)

TIR 70-180 mg/dL	A1C
20%	10.6%
30%	9.8%
40%	9.0%
50%	8.3%
60%	7.5%
70%	6.7%
80%	5.9%
90%	5.1%

**For every 10% increase in TIR
= ~0.8% A1C reduction.**

Beck et al.³

(n = 545 participants with type 1 diabetes)

TIR 70-180 mg/dL	A1C
20%	9.4%
30%	8.9%
40%	8.4%
50%	7.9%
60%	7.4%
70%	7.0%
80%	6.5%
90%	6.0%


**For every 10% increase in TIR
= ~0.5% A1C reduction.**

Post-hoc analysis indicate an association between TIR & diabetes complications



- Post-hoc analysis of SMBG data indicate that each 10-percentage point decrease in TIR was associated with a 64% increase in retinopathy progression and 40% increase in microalbuminuria ($P < 0.001$ for each).¹
- A 10% increase in TIR associated with 19% reduction in urinary albumin excretion.²
- TIR was significantly associated with prevalence of all stages of retinopathy in T2D even after adjusting for clinical risk factors, including A1C.³
- 10% decrease TIR associated with 5% increase in CVD-related mortality and 8% increase in all-cause mortality.⁴

Ambulatory Glucose Profile (AGP) Report and Time in Range (TIR) Inform Therapeutic Decisions

Recommended Time in Range for most people with T1D & T2D¹



DEXCOM | capturAGP™




The y-axis scale and target range are the same as on the AGP graph above. AGP is a summary of glucose values from the report period, with median (50%) and other percentiles shown as if occurring in a single day. Hypothetical patient and AGP report. QR = interquartile range. 1.. Battelino T et al. Diabetes Care. 2019;42(8):1593-1603

Strategies to Help People with Diabetes Increase TIR

Strive for FNIR¹
(flat, narrow, in-range)

Be curious

Discover



Ongoing Learning

**“If _____, then _____.”
Approach**

Use Patterns

Celebrate

Using Strength-Based Language¹ While Talking About TIR

Emotions tied to living with diabetes

Numbers are data, not a judgement*¹

- Not "good" vs "bad" numbers

Develop SMART goals to move towards more TIR and less TBR

Create an action plan together using shared decision-making


- What is one way you might move towards more TIR?
- Would you prefer to eat fewer carbs at breakfast or increase your insulin?

Celebrate the success!

- Identify the positives

G6 Pro report shows a daily TIR bar




3 Joe's best glucose day was April 23, 2021
Joe's glucose data was in the target range about 83% of the day.



* extreme hyper- or hypo-glycemic events warrant serious discussion between a HCP and patient
1. Dickinson, J.K. et al., *Diabetes Care* 2017 Dec; 40(12): 1790-1799.



Frequent* Dexcom CLARITY Views are Associated with up to 15% Increased TIR¹



Frequent* Dexcom CLARITY viewers experience up to 15% increased time spent in range (70-180 mg/dL) as compared to non-users.¹ *Frequent use is defined as four or more monthly log ins to Dexcom CLARITY.
¹ Parker AS, Welsh J, Jimenez A, Walker T. Insights from big data (2): Benefits of self-guided retrospective review of continuous glucose monitoring reports. Diabetes Technol Ther. 2018;20(S1):A-27.

Help Patients Use Dexcom G6 and CLARITY App Settings for More TIR¹⁻³

Dexcom features that may help move towards more **GREEN**

- Wear CGM daily
- Consider incorporating trend arrows into treatment decisions
- Share data with at least one follower¹
- Customize Dexcom G6 high alert settings²
- Use CGM to understand the impact of food, daily activity, stress and medication³


Dexcom features that can reduce **RED**

- Keep Urgent Low Soon Alert enabled to predict when glucose will be less than 55 mg/dL within 20 minutes
- Customize Dexcom G6 low alert settings²

1. Welsh, J et al. Real-Time Sharing and Following of Continuous Glucose Monitoring Data in Youth. Diabetes Ther. <https://doi.org/10.1007/s13300-019-0571-0>

2 Dexcom CGM: Integrating data science and clinical application to support patient glycemic management. ATTD (Advanced Technologies and Treatments in Diabetes) 2020. February 19-22, 2020. Madrid, Spain.

3 Isaacson, B. et al., 2020, Sept 16 *JDST*, Demonstrating the clinical impact of continuous glucose monitoring within an integrated healthcare delivery system. DOI: 10.1177/1932296820955228 e-pub ahead of print



Time in Range in Clinical Practice

Incorporate Discovery Learning with a Shared-Decision Making Approach¹



- Person centered care
- Shared sense of purpose



- Open dialogue based on **clinical evidence and patient preference**



- Encourages active engagement in learning process
- Empowers self-management

- Collaboration between people with diabetes and the care team
- Builds trust



1. https://www.healthit.gov/sites/default/files/nlc_shared_decision_making_fact_sheet.pdf, Accessed March 18, 2021

Dexcom G6 RT-CGM System

Dexcom G6 Overview

The only CGM system **indicated for children aged 2 years and older**

Up to **288 continuous readings per day**

Exceptional accuracy
Class II device designation¹

Zero fingersticks required^{*}

^{*}If your glucose alerts and readings from the G6 do not match symptoms or expectations, use a blood glucose meter to make diabetes treatment decisions



Smart devices sold separately. ^{*}For a list of compatible devices visit www.Dexcom.com/compatibility

Hey Siri, what's my glucose?

Customizable alerts and a fixed Urgent Low alarm
Predictive Urgent Low Soon alert

Data share features with up to **10 followers[†]**

Robust clinical evidence of **improved glycemic outcomes^{2,3}**

Studies prove the clinical benefits of the differentiating attributes of Dexcom.^{4,5}

[†]Separate Follow app required.

1. FDA. <https://www.fda.gov/news-events/press-announcements/fda-authorizes-first-fully-interoperable-continuous-glucose-monitoring-system-streamlines-review>. Accessed June 19, 2019. 2. Beck RW, et al. *JAMA*. 2017;317(4):371-378. 3. Welsh JB, et al. *Diabetes Technol Ther*. 2019;21(3):128-132. 4. Puhr S, et al. *J Diabetes Sci Technol*. 2020;14(1):83-86. 5. Welsh JB, et al. *Diabetes Ther*. 2019;10(2):751-755.


Dexcom G6 Urgent Low Soon Alert Enhances Experience

Smarter Alerts




Urgent Low Soon Alert

- Future-alert function
- Provides earlier actionable alert without increasing nuisance factor



Note: for illustrative purposes only.



CGM = continuous glucose monitoring.

Dexcom CGM: Integrating data science and clinical application to support patient glycemic management. ATTD (Advanced Technologies and Treatments in Diabetes) 2020. February 19-22, 2020. Madrid Spain.

Dexcom


Dexcom CLARITY Allows Patients to Share CGM Data with HCPs and Receive Daily or Weekly Progress Reports on Smartphone¹

Dexcom CLARITY is a cloud-based diabetes management software that helps patients and providers understand and analyze glucose patterns²



Home user²:

- View glucose patterns, trends, and statistics



Healthcare Professional³:

- Graphs show patterns of hypoglycemia and hyperglycemia, allowing providers to prioritize problems and find diabetes management solutions

Dexcom CLARITY Clinic

Simplified Access to Actionable
RT-CGM Data

Get the latest version of CLARITY Uploader. Please install by June 23, 2020.



Welcome to Dexcom CLARITY for Healthcare Professionals

Login with your Dexcom CLARITY Healthcare Professional account

[Forgot your username?](#) [Forgot your password?](#)

Need to register your clinic? [Register Now](#)

View Data from a Dexcom Receiver as Guest

Upload a Dexcom CGM receiver without logging in. This one-time upload will allow you to view a report for the data from your CGM receiver only, but the data will not be saved to your account.

View data shared from a smart device

If your patient has the [Dexcom CLARITY app](#) on their smart device, they can generate a data-sharing code so you can view their data on your schedule.

Enter patient provided sharing code

Not the best way to set up your CLARITY clinic account!

Patient list- add new patients in < 1 minute, but do need DOB

Search for Patient name or Patient ID Add new patient

E	DOB	PATIENT ID	LAST UPLOADED	DATA SHARING
	Feb 2, 1947		Jun 21, 2020	✓ On
	Nov 17, 1936		Mar 30, 2020	✓ On
	Apr 29, 1946		Mar 30, 2020	✓ On
	Dec 22, 1960		Dec 21, 2018	✗ Off
	Apr 15, 1966		Jun 21, 2020	✓ On
	Oct 28, 1999		Jun 21, 2020	✓ On
	Jul 16, 1986		Feb 25, 2020	✓ On
	Jul 16, 1986		Aug 9, 2018	✗ Off
	Sep 7, 1957		Jun 21, 2020	✓ On

Our Patient Named “Dexcom Clarity,” DOB 1-1-2015: Set up an Indefinite Share Quickly

Get the latest version of CLARITY Uploader. Please install by June 23, 2020.

Dexcom CLARITY Patients Settings Support Nicholas Argento, MD
Maryland Endocrine ...

clarity ⊗ Add new patient Export all data

PATIENT NAME	DOB	PATIENT ID	LAST UPLOADED	DATA SHARING
Clarity, Dexcom	Jan 1, 2015			✖ Off ⊗


Upload data Save or print report Go to interactive reports

Delete Edit Export Share data

Dexcom © 2016-2020 Contact Us Terms of Use Privacy Policy Patient Bill of Rights Safety Information

Dexcom and Dexcom CLARITY are registered trademarks of Dexcom, Inc. in the United States and may be in other countries. All other marks are property of their respective owners.

00386270000491 • Dexcom CLARITY v3.28.1 • PN 350-0011 • DOM 2020-05-26



Invite Patient to Share Data

The screenshot shows a web browser window displaying the Dexcom CLARITY Clinic Portal. The URL is clarity.dexcom.com/professional/patients/1670144165472677888/share. The page features a navigation bar with 'Patients', 'Settings', and 'Support' options. A user profile for 'Nicholas Argento, MD, Maryland Endocrine ...' is visible in the top right. A central dialog box titled 'Invite this patient to share data' is open, containing the following text: 'If the patient accepts, their personal Dexcom CLARITY account and your clinic's account will automatically share data between them.' Below this, a patient profile for 'Dexcom Clarity' is shown with a 'DOB' of 'Jan 1, 2015' and a 'Patient ID' field. The dialog prompts the user to 'Please select one of the following options.' with radio buttons for 'Print an Invitation' (selected) and 'Email an Invitation'. At the bottom of the dialog are 'Invite' and 'Cancel' buttons. The background shows a patient record with fields for 'PATIENT NAME' (Clarity, Dexcom) and buttons for 'Export all data' and 'Share data'. The footer of the page includes copyright information: 'Dexcom and Dexcom CLARITY are registered trademarks of Dexcom, Inc. in the United States and may be in other countries. All other marks are property of their respective owners.' and version information: '00386270000491 • Dexcom CLARITY v3.28.1 • PN 350-0011 • DOM 2020-05-26'.

- After clicking Share data option, click invite, print or email
- Patient needs to have the Dexcom CLARITY app on their phone or computer to enter code.


Share Data with Clinic

The screenshot shows a web browser window with the URL `clarity.dexcom.com/professional/patients/1670144165472677888/share/print?autoPrint=true`. The page displays a sharing code generation interface for a patient named "Dexcom CLARITY". The code generated is **HSZL-XPSC-CFVT**. The interface includes the following text:

- Patient Name:** Dexcom CLARITY
- Generated at:** Jun 21, 2020 5:08 PM
- Dexcom CLARITY** logo
- Share data with your clinic**
- Maryland Endocrine PA invites you to share your data using Dexcom CLARITY.**
- Dexcom CLARITY software captures your continuous glucose monitoring (CGM) data so you and your clinic can view patterns, trends and statistics anytime, anywhere.
- Your sharing code***
- HSZL-XPSC-CFVT**
- Enter this code at <https://clarity.dexcom.com/share>
- *Expires: July 21, 2020
- See all your CGM data on the go with the Dexcom CLARITY app. Visit <https://clarity.dexcom.com/mobile> to get started.
- Buttons: Print, Close

- Ask patient to enter code in their Dexcom CLARITY account under Profile
- Then authorize sharing, accept code
- (not generate code)

Case Study



Case Study #1

53 yo man with long standing diabetes presents as a new patient to DM clinic. He does not check his FSBG on a regular basis.

- MF 2000 mg qd
- Tresiba 10 units qd
- HbA1c 13.2 %

No changes to medication, placed a Dexcom G6 during office apt.

Glucose Statistics


Avg Glucose mg/dL

243


Glucose Exposure

Very Low	Low	In Target Range	High	Very High
< 54 mg/dL	< 70 mg/dL	70 - 180 mg/dL	> 180 mg/dL	> 250 mg/dL
0.0%	0.0%	15.0%	85.0%	39.3%
Glucose Ranges				

Coefficient of Variation	SD mg/dL
25.1%	61
Glucose Variability	



Curves/plots represent glucose frequency distributions by time regardless of date.





Case Study #1


At follow up apt after 2 weeks meal time insulin Humalog started.

- Humalog 10 units TID with a 1:50 correction about 150 mg/dL


Patient returns one month later for a follow up apt.

Compare Report

Baseline



One Month Later



Polling question 4

- Which patient do you think would be the best candidate for RT-CGM?
 - A. 30 y.o. man with type 1 diabetes and A1C 7.5%
 - B. 55 y.o. woman with type 2 diabetes, A1C 9% on basal insulin
 - C. 65 y.o. man with type 1 diabetes with frequent hypoglycemia
 - D. 40 y.o. woman with type 2 diabetes, A1C 9% on metformin and Sulfonylurea
 - E. All of the above

Summary

TAKEAWAY

- Metrics beyond A1C are needed to fully understand the variation in glucose values

GOAL:

- The goal is to increase time in range (TIR) while decreasing time below range

ACRONYMS TO REMEMBER:

- **MGLR** or more green, less red on the TIR bar
- **FNIR** or flat, narrow, in range on the trends graph

ADDITIONAL

- TIR discussions can focus on strengths and making positive changes for health
- Avid Dexcom CLARITY users* on average, spend 15% more TIR (70-180 md/dL that non-users)¹

*Avid users have 4 or more monthly long-ins to Dexcom CLARITY. 1 Parker AS, Welsh J, Jimenez A, Walker T. Insights from big data (2): Benefits of self-guided retrospective review of continuous glucose monitoring reports. Diabetes Technol Ther. 2018;20(S1):A-27.

Polling question 5

- How likely are you to discuss RT-CGM with a patient tomorrow?
 - A. Unlikely, I need more information
 - B. I think I will, I'm unsure about the time it might take
 - C. I definitely will, I believe in the value of RT-CGM
 - D. I'm already discussing RT-CGM every day

Resources to Download

Why Time in Range Matters?


dexcom

WHAT IS TIME IN RANGE (TIR)?

Time in Range or TIR is the percentage of time glucose levels are in target range. For most people with diabetes the target range is between 70 and 180 mg/dL.†

The daily goal for most people with type 1 diabetes (T1D) and type 2 diabetes (T2D) is >70% TIR or approximately 16 hours, 48 min.†

See the TIR bar for recommended targets.



RECOMMENDED TIME IN RANGE

Recommended Time in Range for most people with T1D & T2D†



Goal: Increase TIR while decreasing TBR
See **More Green** and **Less Red (MGLR)** on the TIR bar



Each 5% increase in TIR is considered clinically significant.¹
For every 10% increase in TIR = ~0.8% A1C reduction.²

Metrics Beyond A1C

dexcom

THE VALUE OF METRICS BEYOND A1C


A1C is the average glucose over the last 2-3 months but does not identify glycemic variability¹

CGM data can identify patterns of hypo- and hyperglycemia, assess glycemic excursions and glucose variability to allow for therapy modification.¹



Same A1C but CGM Patterns Drive Different Treatment Plans²

Estimated A1C for a Time in Range (TIR) level³



Each 5% increase in TIR is considered clinically significant.¹
For every 10% increase in TIR = ~0.8% A1C reduction.³



KEY METRICS

Number of days with CGM data
14+ days recommended


Percentage of time CGM is active
>70% of data recommended

Mean glucose
The average glucose

Glucose Management Indicator (GMI)
Approximate A1C levels based on average glucose measured using CGM values

Coefficient of Variation (CV)
Measure of glycemic variability ≤36% is recommended¹

GOALS FOR TIR



dexcom

OUR MISSION, TOGETHER

Too many people with diabetes are still unaware of modern diabetes metrics, such as time in range.

We want to change that.

#WHENINRANGE



dexcom

Beyond Type 1

JDRF IMPROVING LIVES. CURING TYPE 1 DIABETES.

Children with Diabetes®
www.childrenwithdiabetes.com

TCOYD
TAKING CONTROL OF YOUR DIABETES



COLLEGE DIABETES NETWORK

Dexcom is Here to Support YOU and Your Patients



Customer Sales Support

Support with Dexcom orders and general customer questions

1-888-738-3646

[Place Your First Order](#)

[Need Pharmacy Information](#)



Global Technical Support

Product troubleshooting or replacement inquiries

1-844-607-8398

Available 24 hours a day; 7 days a week

[Request a Call Back](#)

[Submit a Patient Support Request](#)

[Request Sensor Overpatches](#)

[Chat Live with Dexcom Tech Support](#)



Dexcom CARE

Dexcom CGM training, software downloads, and tutorials

1-877-339-2664



See dexcom.com/contact for current contact hours

Connect with Global Thought Leaders to Advance Diabetes Technology



Visit SIGMA to learn about:

- Expert Presentations
- Clinical Evidence
- Diabetes Technology Resources
- Coverage and Reimbursement
- Latest News



To request membership,
scan the QR code

If you are not already a member, we invite you to learn more about SIGMA or request membership at www.cgmonitoring.net

Dexcom G6 Safety Statement

Failure to use the Dexcom G6 Continuous Glucose Monitoring System (G6) and its components according to the instructions for use provided with your device and available at <https://www.dexcom.com/safety-information> and to properly consider all indications, contraindications, warnings, precautions, and cautions in those instructions for use may result in you missing a severe hypoglycemia (low blood glucose) or hyperglycemia (high blood glucose) occurrence and/or making a treatment decision that may result in injury. If your glucose alerts and readings from the G6 do not match symptoms or expectations or you're taking over the recommended maximum dosage amount of 1000mg of acetaminophen every 6 hours, use a blood glucose meter to make diabetes treatment decisions. Seek medical advice and attention when appropriate, including for any medical emergency.

The web-based Dexcom CLARITY software is intended for use by both home users and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis, and evaluation of historical CGM data to support effective diabetes management. It is intended for use as an accessory to Dexcom CGM devices with data interface capabilities. Caution: The software does not provide any medical advice and should not be used for that purpose. Home users must consult a healthcare professional before making any medical interpretation and therapy adjustments from the information in the software. Caution: Healthcare professionals should use information in the software in conjunction with other clinical information available to them. Caution: Federal (US) law restricts this device to sale by or on the order of a licensed healthcare professional.

Dexcom, Dexcom G6, Dexcom Share, and Dexcom CLARITY are registered trademarks of Dexcom, Inc. in the United States and/or other countries.

Dexcom, Inc. | 6340 Sequence Drive | San Diego, CA 92121

Technical Support: 1-877-339-2664 | www.dexcom.com

©2021 Dexcom, Inc. All rights reserved.

Dexcom G6 Pro Safety Statement

Failure to use the Dexcom G6 Pro Continuous Glucose Monitoring System (G6 Pro) and its components according to the instructions for use provided with your device and available at <https://www.dexcom.com/safety-information> and to properly consider all indications, contraindications, warnings, precautions, and cautions in those instructions for use may result in your patient missing a severe hypoglycemia (low blood glucose) or hyperglycemia (high blood glucose) occurrence and/or making a treatment decision that may result in injury. If your patient's glucose alerts and readings from the G6 Pro do not match symptoms or expectations or your patient is taking over the recommended maximum dosage amount of 1000mg of acetaminophen every 6 hours, use a blood glucose meter to make diabetes treatment decisions. Your patient will not receive alerts and alarms when the G6 Pro is on blinded mode. Seek medical advice and attention when appropriate, including for any medical emergency.

The web-based Dexcom CLARITY software is intended for use by both home users and healthcare professionals to assist people with diabetes and their healthcare professionals in the review, analysis, and evaluation of historical CGM data to support effective diabetes management. It is intended for use as an accessory to Dexcom CGM devices with data interface capabilities. Caution: The software does not provide any medical advice and should not be used for that purpose. Home users must consult a healthcare professional before making any medical interpretation and therapy adjustments from the information in the software. Caution: Healthcare professionals should use information in the software in conjunction with other clinical information available to them. Caution: Federal (US) law restricts this device to sale by or on the order of a licensed healthcare professional.

Dexcom, Dexcom G6, Dexcom Share, and Dexcom CLARITY are registered trademarks of Dexcom, Inc. in the United States and/or other countries.

Dexcom, Inc. | 6340 Sequence Drive | San Diego, CA 92121
Technical Support: 1-877-339-2664 | www.dexcom.com

THANK YOU!!

Questions?



Smart devices sold separately.